

Salt Water



802 Washington Ave. • Chestertown • Maryland • 21620 • USA
800-344-3100 • +1 410-778-3100 • Fax 410-778-6394
Visit us on the web at



WARNING: Cancer and Reproductive Harm -
www.P65Warnings.ca.gov

Test Procedure

1. Fill the tube [0715] to the 5 mL line with the sample water.
2. Add three drops of *Chloride Reagent #1 [4504] to sample. Swirl to mix. Solution will turn yellow.
3. Holding *Salt Titrant [3741LWT] in a vertical position, add one drop at a time, swirling to mix, until the color changes from yellow to milky orange. Count the number of drops required to reach the final color. Add titrant slowly, allowing each drop to form fully on the tip before dispensing. Rapid titration with incomplete drops will result in erroneous high readings.
4. Multiply number of drops used by 200. Record result as ppm salt [sodium chloride].

*Salt Titrant will stain. Avoid contact with skin. Place tip of titrant bottle well into titration tube while dispensing.



Store the kit at room temperature.

*Potential health hazard. Read SDS at lamotte.com.

Emergency information for all LaMotte reagents is available from Chem-Tel [US: 1-800-255-3924]
[International: call collect, 813-248-0585].

© 2022 LaMotte Company, All Rights Reserved | 63156 | 02.11.2022

