

VEE GEE

GLASSWARE & ACCESSORIES

VEE GEE Scientific
Laboratory Equipment & Supplies



1.888.610.7664

 www.calcert.com

sales@calcert.com

Separatory Funnels

Feature smooth tooled tips, large white block letters, and large marking squares

- Supplied with PTFE stopcocks to eliminate freezing and lubricant contamination
- Made from high-quality borosilicate tubing with a low coefficient of expansion

SIBATA Separatory Funnels with Stopper - Squibb Style

SIBATA SEPARATORY FUNNELS WITH STOPPERS - SQUIBB STYLE								
Cat. No.	Capacity (mL)	Stopper No.	Stopcock Plug Size	Approx. O.D. x Stem Length (in)	Approx. O.D. x Stem Length (mm)	Approx. Total Length (in/mm)	Pack Qty.	Case Qty.
1426T-300	300	24/40	4	0.4 x 2.8	10 x 70	15.4 / 390	1	6
1426T-500	500	24/40	4	0.4 x 2.8	10 x 70	17 / 430	1	6
1426T-1000	1000	29/42	6	0.5 x 3.1	12 x 80	19.3 / 490	1	4
1426T-2000	2000	29/42	6	0.5 x 3.1	12 x 80	20 / 510	1	2



VEE GEE Separatory Funnels with Stopper - Squibb Style

- All sizes supplied with a polyethylene standard taper stopper
- Meet ASTM E-1096 Type 4 Specifications for Classification, Design, and Markings

VEE GEE SEPARATORY FUNNELS WITH STOPPERS - SQUIBB STYLE								
Cat. No.	Capacity (mL)	Stopper No.	Stopcock Plug Size	Approx. O.D. x Stem Length (in)	Approx. O.D. x Stem Length (mm)	Approx. Total Length (in/mm)	Pack Qty.	Case Qty.
20149-250	250	22	2	0.4 x 2.6	10 x 65	13.6 / 345	1	24
20149-500	500	27	4	0.4 x 2.6	10 x 65	15.4 / 392	1	24

