

VEE GEE

GLASSWARE & ACCESSORIES

VEE GEE Scientific
Laboratory Equipment & Supplies



1.888.610.7664

 www.calcert.com

sales@calcert.com

Lead Rings

Stabilize vessels and labware on the bench or in water baths

- Vinyl-coated lead rings will not mark or scratch benchtops and labware
- C-shape format allows some manual size variation for a customized fit

LEAD RINGS								
Cat No.	Shape	Color	Flask Sizes (mL)	ID (in)	ID (cm)	OD (in)	OD (cm)	Weight* (g)
20335-1	C	Blue	125–500	1.9	4.2	3	7.5	215
20335-2	O	Blue	125–500	1.9	4.8	3	7.7	210
20335-3	C	Blue	250–1000	2	5.1	3.5	9	465
20335-4	O	Blue	250–1000	2	5.1	3.5	9	490
20335-5	C	Blue	500–2000	2.1	5.4	3.9	9.8	595
20335-6	O	Blue	500–2000	2.2	5.7	3.9	9.8	610
20335-7	C	Blue	1000–4000	2.6	6.6	4.4	11.3	890
20335-8	O	Blue	1000–4000	2.8	7	4.5	11.5	935

**Weights may vary by +/- 3%*



Preserve your solutions with
Parafilm® M self-sealing film *pg. 14*

