

## Oakton® Electrode Cleaning and Storage Solutions

### Clean electrodes make clean measurements

- Easily remove grease, protein, and other deposits with cleaning solution
- Make electrode last longer—store in our storage solution
- Choose 500-mL bottle or case of 24 500-mL bottles



**OAKTON®**

| Solution type | Catalog number              | Price/ea | Catalog number              | Price/cs |
|---------------|-----------------------------|----------|-----------------------------|----------|
| Cleaning      | <a href="#">GY-00653-06</a> |          | <a href="#">GY-35655-06</a> |          |
| Storage       | <a href="#">GY-00653-04</a> |          | <a href="#">GY-35655-04</a> |          |

## Oakton® Reference Fill Solutions

### Effortless refilling!

- Choose your solution: 3.8 KCl with AgCl for Ag/AgCl electrodes or 4 KCl for calomel
- Sold in packs of two 125-mL bottles



**OAKTON®**

| Description       | Catalog number              | Price/pk |
|-------------------|-----------------------------|----------|
| 3.8 KCl with AgCl | <a href="#">GY-05992-47</a> |          |
| 4 KCl             | <a href="#">GY-05992-48</a> |          |

## Electrode Storage Solution

### Optimize your ROSS™ electrode performance

- Specially formulated for ROSS pH electrodes
- 475-mL bottle



| Description                        | Catalog number              | Price |
|------------------------------------|-----------------------------|-------|
| ROSS pH electrode storage solution | <a href="#">GY-05730-60</a> |       |

## Oakton® Grip-Clip™ Holder

### Hands-free measurements

**NEW**

- No messing with electrode holders—simply clip to side of beaker for easier probe handling
- Holds an ATC probe and a 11- or 12-mm probe or a 16-mm probe



**OAKTON®**

| Description             | Catalog number              | Price |
|-------------------------|-----------------------------|-------|
| Oakton Grip-Clip holder | <a href="#">GY-35614-50</a> |       |

## pH Electrode Cleaning Kit

### Get longer electrode life

- Restores performance of your electrode
- Remove protein, bacteria, oil, and grease with ease
- Includes one bottle each of cleaning solution, beaker, and pipette



**NEW**

| Description               | Catalog number              | Price |
|---------------------------|-----------------------------|-------|
| pH electrode cleaning kit | <a href="#">GY-05730-70</a> |       |

#### Accessories

**GY-05730-71 Replacement solution A** for removal of proteins. Four 30-mL bottles, beaker, and pipette

**GY-05730-72 Replacement solution B** for removal of bacteria. Four 60-mL bottles, beaker, and pipette

**GY-05730-73 Replacement solution C** for general-purpose cleaning. Four 30-mL bottles, beaker, and pipette

**GY-05730-74 Replacement solution D** for removal of oil and grease. Four 60-mL bottles, beaker, and pipette

## Electrode Saver Bottles

### Keeps electrodes in working order

- Protect against bumps, spills, and drying out
- Accepts electrodes up to 13 mm diameter

| Capacity (mL) | Dimensions (H x dia)      | Catalog number              | Price |
|---------------|---------------------------|-----------------------------|-------|
| 50            | 4" x 1½" (10 x 3.8 cm)    | <a href="#">GY-05990-99</a> |       |
| 5             | 1¾" x 7/8" (4.4 x 2.2 cm) | <a href="#">GY-35805-50</a> |       |



05990-99

35805-50

## Electrode Storage Sleeve

### Eliminate broken pH bulbs

- Protect fragile pH glass bulbs up to 12-mm dia
- Keep electrodes conditioned and ready for pH measurement
- Weighted base for safe, upright storage

**NEW**



| Description              | Catalog number              | Price |
|--------------------------|-----------------------------|-------|
| Electrode storage sleeve | <a href="#">GY-99570-50</a> |       |