



Bio Plas G-Tube®

Premium Flat Top Microcentrifuge Tubes

The **G-Tube™** Flat Top Microcentrifuge Tube is virtually transparent and far exceeds the toughness of any microcentrifuge tubes. When the **G-Tube™** Flat Top Microcentrifuge Tube cap is closed there will be an audible "snap" that indicates a secure closure has been formed.

Made of tough, clear, proprietary formulated polypropylene, the 0.6mL, 1.5mL and 2.0mL

G-Tube™ has the same chemical resistant quality as standard tubes, but they're sturdier and virtually transparent to give a better view. The flat top cap is easy to open and still provides a maximum secure seal. All three volumes are available in siliconized material. Siliconized **G-Tubes™** are not recommended for autoclaving.

The **G-Tube™** is built for high impact usage and will withstand over 30,000 x g when spun. Thicker walls make this the tube of choice for extreme stress, high spinning speeds, autoclaving, boiling, chemical resistance and freezing. Will withstand temperature ranges of -80°C (-112°F) to 121 °C (250°F).



0.6mL **G-Tube™** Flat Top Microcentrifuge Tube, Green, 4044



BIO PLAS, INC. A WOMAN OWNED SMALL BUSINESS

4380 Redwood Hwy, Suite C-16, San Rafael, CA 94903 USA

Tel: (415) 472-3777

Fax: (415) 472-3758

info@bioplas.com

www.bioplas.com



BIO PLAS INC.
 INNOVATIVE IDEAS...
 FOR TODAY'S LABORATORY DISPOSABLES
 MADE IN THE USA



0.6mL G-Tube® Flat Top
Microcentrifuge Tube, Yellow, 4042



1.5mL G-Tube® Flat Top
Microcentrifuge Tube, Blue 4031



2.0mL G-Tube® Flat Top
Microcentrifuge Tube, Green, 4056

G-Tube™ Flat Top Microcentrifuge Tube Features:

- Made from high quality medical grade resins
- RNase and DNase free. Human DNA free. Non-Pyrogenic
- Ultra low retention
- Lid Closure: The flat top lid forms an extra-tight seal that will not come open during processing. Yet, it's easy to open with thumb pressure or our handy MicroTube Tool (catalog #4000.)
- Bio Plas G-Tubes have a textured writing surface on top of the caps. The 1.5ml and 2.0 mL tubes have an additional textured writing surface on the side of the tubes
- Graduation: The 1.5ml tubes and 2.0ml tubes are graduated every 0.1ml
- Boilproof
- Will withstand speeds up to 30000 X G
- Easy needle access through the piercing port in the exact center of the cap
- Autoclavable at 250°F at 15 psi for 15 minutes.

Ordering Information					
0.6 mL Catalog #	1.5 mL Catalog #	2.0 mL Catalog #	Description	Pkg. Qty	Pkgs./Case
4040	4030	4050	G-Tube® Flat Top Microcentrifuge Tube, Natural	500	20
4041	4031	4051	G-Tube® Flat Top Microcentrifuge Tube, Red	500	20
4042	4032	4052	G-Tube® Flat Top Microcentrifuge Tube, Yellow	500	20
4043	4033	4053	G-Tube® Flat Top Microcentrifuge Tube, Blue	500	20
4044	4034	4054	G-Tube® Flat Top Microcentrifuge Tube, Green	500	20
4046	4036	4056	G-Tube® Flat Top Microcentrifuge Tube, Orange	500	20
4047	4037	4057	G-Tube® Flat Top Microcentrifuge Tube, Lavender	500	20
4048	4038	4058	G-Tube® Flat Top Microcentrifuge Tube, Amber	500	20
4039	4049	4059	G-Tube® Flat Top Microcentrifuge Tube, Rainbow	500	20

BIO PLAS, INC. A WOMAN OWNED SMALL BUSINESS

4380 Redwood Hwy, Suite C- 116, San Rafael, CA 94903 USA Tel: (415) 472-3777 Fax: (415) 472-3758

info@bioplas.com www.bioplas.com

G-Tube® is a registered trademark of Bio Plas, Inc. San Rafael, CA

1.888.610.7664



www.calcert.com

sales@calcert.com